Package 'affyPLM'

May 11, 2015

2 3 6 7 8 9 11.

bg.correct.LESN

bg.c	orrect.LESN	LE	SN	- L	юи	, E	nd	Si	gn	ıai	l i.	s N	Vo	ise	B	ac	cks	gra	ои	nd	cc	ori	ec	cti	on	s			
Index																													21
	threestepPLM							•																					19
	threestep																												
	ReadRMAExpress rmaPLM																												
	pseudo.coloring .																												

Description

This function background corrects PM probe data using LESN - Low End Signal is Noise concepts.

Usage

```
bg.correct.LESN(object, method=2, baseline=0.25, theta=4)
```

Arguments

object an AffyBatch

method an integer code specifying which method to use

baseline A baseline value to use

theta A parameter used in the background correction process

Details

This method will be more formally documented at a later date.

The basic concept is to consider that the lowest end of intensites is most likely just noise (and should be heavily corrected) and the highest end signals are most likely signal and should have little adjustment. Low end signals are made much smaller while high end signals get less adjustment relative adjustment.

Value

An AffyBatch

Author(s)

Ben Bolstad <bmb@bmbolstad.com>

References

Bolstad, BM (2004) Low Level Analysis of High-density Oligonucleotide Array Data: Background, Normalization and Summarization. PhD Dissertation. University of California, Berkeley.

fitPLM 3

Examples

```
if (require(affydata)) {
  data(Dilution)
  Dilution.example.bgcorrect <- bg.correct.LESN(Dilution)
}</pre>
```

fitPLM

Fit a Probe Level Model to Affymetrix Genechip Data.

Description

This function converts an AffyBatch into an PLMset by fitting a specified robust linear model to the probe level data.

Usage

Arguments

object an AffyBatch

model A formula describing the model to fit. This is slightly different from the standard

method of specifying formulae in R. Read the description below

variable.type a way to specify whether variables in the model are factors or standard variables

constraint.type

should factor variables sum to zero or have first variable set to zero (endpoint

constraint)

subset a vector with the names of probesets to be used. If NULL then all probesets are

used.

normalize logical value. If TRUE normalize data using quantile normalization

background logical value. If TRUE background correct using RMA background correction

background.method

name of background method to use.

normalize.method

name of normalization method to use.

background.param

A list of parameters for background routines

4 fitPLM

```
normalize.param
```

A list of parameters for normalization routines

output.param A list of parameters controlling optional output from the routine.

model.param A list of parameters controlling model procedure

verbosity.level

An integer specifying how much to print out. Higher values indicate more verbose. A value of 0 will print nothing

Details

This function fits robust Probe Level linear Models to all the probesets in an AffyBatch. This is carried out on a probeset by probeset basis. The user has quite a lot of control over which model is used and what outputs are stored. For more details please read the vignette.

Value

An PLMset

Author(s)

Ben Bolstad <bmb@bmbolstad.com>

References

Bolstad, BM (2004) Low Level Analysis of High-density Oligonucleotide Array Data: Background, Normalization and Summarization. PhD Dissertation. University of California, Berkeley.

See Also

```
expresso, rma, threestep
```

Examples

```
if (require(affydata)) {
   data(Dilution)
   Pset <- fitPLM(Dilution, model=PM ~ -1 + probes + samples)
   se(Pset)[1:5,]

image(Pset)
   NUSE(Pset)

#now lets try a wider class of models
   ## Not run: Pset <- fitPLM(Dilution, model=PM ~ -1 + probes + liver,
   normalize=FALSE, background=FALSE)

## End(Not run)
   ## Not run: Coefs(Pset)[1:10,]

## Not run: Pset <- fitPLM(Dilution, model=PM ~ -1 + probes + liver +
   scanner, normalize=FALSE, background=FALSE)

## End(Not run)
   coefs(Pset)[1:10,]</pre>
```

MAplot 5

```
#try liver as a covariate
logliver <- log2(c(20,20,10,10))
## Not run: Pset <- fitPLM(Dilution, model=PM~-1+probes+logliver+scanner,
normalize=FALSE, background=FALSE, variable.type=c(logliver="covariate"))
## End(Not run)
coefs(Pset)[1:10,]

#try a different se.type
## Not run: Pset <- fitPLM(Dilution, model=PM~-1+probes+scanner,
normalize=FALSE,background=FALSE,m odel.param=list(se.type=2))
## End(Not run)
se(Pset)[1:10,]
}</pre>
```

MAplot

Relative M vs. A plots

Description

Create boxplots of M or M vs A plots. Where M is determined relative to a specified chip or to a pseudo-median reference chip.

Arguments

A Additional parameters for the routine

A A vector to plot along the horizonal axis

M A vector to plot along vertical axis

subset A set of indices to use when drawing the loess curve show.statistics

If true some summary statistics of the M values are drawn span span to be used for loess fit.

family.loess "guassian" or "symmetric" as in loess.

cex Size of text when writing summary statistics on plot

See Also

```
mva.pairs
```

```
normalize.ExpressionSet
```

Normalization applied to ExpressionSets

Description

Allows the user to apply normalization routines to ExpressionSets.

Usage

Arguments

eset	An ExpressionSet
span	parameter to be passed to the function loess.
choose.subset	use a subset of values to establish the normalization relationship
subset.size	number to use for subset
verbose	verbosity flag
family	parameter to be passed to the function loess.
prd.td	cutoff parameter (details in the bibliographic reference)
trim	How much to trim from the top and bottom before computing the mean when using the scaling normalization
baseline	Index of array to use as baseline, negative values (-1,-2,-3,-4) control different baseline selection methods
transfn	Transform the ExpressionSet before normalizing. Useful when dealing with expression values that are log-scale
baseline.type	A method of selecting the baseline array
	Additional parameters that may be passed to the normalization routine

Details

This function carries out normalization of expression values. In general you should either normalize at the probe level or at the expression value level, not both.

Typing normalize. ExpressionSet.methods should give you a list of methods that you may use. note that you can also use the normalize function on ExpressionSets. Use method to select the normalization method.

Value

A normalized ExpressionSet.

Author(s)

Ben Bolstad, <bmb@bmbolstad.com>

References

Bolstad, BM (2004) Low Level Analysis of High-density Oligonucleotide Array Data: Background, Normalization and Summarization. PhD Dissertation. University of California, Berkeley.

Examples

```
if (require(affydata)) {
  data(Dilution)
  eset <- rma(Dilution, normalize=FALSE, background=FALSE)
  normalize(eset)
}</pre>
```

normalize.quantiles.probeset

Quantile Normalization applied to probesets

Description

Using a normalization based upon quantiles, this function normalizes a matrix of probe level intensities.

Usage

```
normalize.AffyBatch.quantiles.probeset(abatch,type=c("separate","pmonly","mmonly","together"),use
```

Arguments

```
abatch An AffyBatch
```

type how should MM and PM values be handled

use.median use median rather than mean use.log take logarithms, then normalize

8 normalize.scaling

Details

This function applies the quantile method in a probeset specific manner.

In particular a probeset summary is normalized using the quantile method and then the probes adjusted accordingly.

Value

A normalized AffyBatch.

Author(s)

Ben Bolstad, <bmb@bmbolstad.com>

References

Bolstad, B (2001) *Probe Level Quantile Normalization of High Density Oligonucleotide Array Data.* Unpublished manuscript http://oz.berkeley.edu/~bolstad/stuff/qnorm.pdf

Bolstad, B. M., Irizarry R. A., Astrand, M, and Speed, T. P. (2003) *A Comparison of Normalization Methods for High Density Oligonucleotide Array Data Based on Bias and Variance*. Bioinformatics 19(2) ,pp 185-193. http://www.stat.berkeley.edu/~bolstad/normalize/normalize.html

See Also

```
normalize.quantiles
```

normalize.scaling

Scaling normalization

Description

Allows the user to apply scaling normalization.

Usage

Arguments

X A matrix. The columns of which are to be normalized.

abatch An AffyBatch

type A parameter controlling how normalization is applied to the Affybatch.

trim How much to trim from the top and bottom before computing the mean when

using the scaling normalization.

PLMset-class 9

baseline Index of array to use as baseline, negative values (-1,-2,-3,-4) control different baseline selection methods.

log.scalefactors

Compute the scale factors based on log2 transformed data.

Details

These function carries out scaling normalization of expression values.

Value

A normalized ExpressionSet.

Author(s)

Ben Bolstad, <bmb@bmbolstad.com>

Examples

```
if (require(affydata)) {
  data(Dilution)
  normalize.AffyBatch.scaling(Dilution)
}
```

PLMset-class

Class PLMset

Description

This is a class representation for Probe level Linear Models fitted to Affymetrix GeneChip probe level data.

Objects from the Class

Objects can be created using the function fitPLM

Slots

probe.coefs: Object of class "matrix". Contains model coefficients related to probe effects.

se.probe.coefs: Object of class "matrix". Contains standard error estimates for the probe coefficients.

chip.coefs: Object of class "matrix". Contains model coefficients related to chip (or chip level) effects for each fit.

se.chip.coefs: Object of class "matrix". Contains standard error estimates for the chip coefficients.

const.coefs: Object of class "matrix". Contains model coefficients related to intercept effects for each fit.

10 PLMset-class

se.const.coefs: Object of class "matrix". Contains standard error estimates for the intercept estimates

model.description: Object of class "character". This string describes the probe level model fitted.

weights: List of objects of class "matrix". Contains probe weights for each fit. The matrix has columns for chips and rows are probes.

phenoData: Object of class "phenoData" This is an instance of class phenoData containing the patient (or case) level data. The columns of the pData slot of this entity represent variables and the rows represent patients or cases.

annotation A character string identifying the annotation that may be used for the ExpressionSet instance.

experimentData: Object of class "MIAME". For compatibility with previous version of this class description can also be a "character". The class characterOrMIAME has been defined just for this

cdfName: A character string giving the name of the cdfFile.

nrow: Object of class "numeric". Number of rows in chip.

ncol: Object of class "numeric". Number of cols in chip.

narrays: Object of class "numeric". Number of arrays used in model fit.

normVec: Object of class "matrix". For storing normalization vector(s). Not currentl used

varcov: Object of class "list". A list of variance/covariance matrices.

residualSE: Object of class "matrix". Contains residual standard error and df.

residuals: List of objects of class "matrix". Contains residuals from model fit (if stored).

model.call: Object of class "call"

Methods

```
weights<- signature(object = "PLMset"): replaces the weights.</pre>
```

weights signature(object = "PLMset"): extracts the model fit weights.

coefs<- signature(object = "PLMset"): replaces the chip coefs.</pre>

coefs signature(object = "PLMset"): extracts the chip coefs.

se signature(object = "PLMset"): extracts the standard error estimates of the chip coefs.

se<- signature(object = "PLMset"): replaces the standard error estimates of the chip coefs.

coefs.probe signature(object = "PLMset"): extracts the probe coefs.

se.probe signature(object = "PLMset"): extracts the standard error estimates of the probe coefs.

coefs.const signature(object = "PLMset"): extracts the intercept coefs.

se.const signature(object = "PLMset"): extracts the standard error estimates of the intercept coefs.

getCdfInfo signature(object = "PLMset"): retrieve the environment that defines the location of probes by probe set.

PLMset2exprSet 11

image signature(x = "PLMset"): creates an image of the robust linear model fit weights for each sample.

indexProbes signature(object = "PLMset", which = "character"): returns a list with
 locations of the probes in each probe set. The list names defines the probe set names. which
 can be "pm", "mm", or "both". If "both" then perfect match locations are given followed by
 mismatch locations.

Mbox signature(object = "PLMset"): gives a boxplot of M's for each chip. The M's are computed relative to a "median" chip.

normvec signature(x = "PLMset"): will return the normalization vector (if it has been stored).

residSE signature(x = "PLMset"): will return the residual SE (if it has been stored).

boxplot signature(x = "PLMset"): Boxplot of Normalized Unscaled Standard Errors (NUSE).

NUSE signature(x = "PLMset"): Boxplot of Normalized Unscaled Standard Errors (NUSE) or NUSE values.

RLEI signature(x = "PLMset"): Relative Log Expression boxplot or values.

Note

This class is better described in the vignette.

Author(s)

References

Bolstad, BM (2004) Low Level Analysis of High-density Oligonucleotide Array Data: Background, Normalization and Summarization. PhD Dissertation. University of California, Berkeley.

PLMset2exprSet

Convert a PLMset to an ExpressionSet

Description

This function converts a PLMset to an ExpressionSet. This is often useful since many Bioconductor functions operate on ExpressionSet objects.

Usage

```
PLMset2exprSet(pset)
pset2eset(pset)
```

Arguments

pset

The PLMset to convert to ExpressionSet.

12 preprocess

Details

These functions convert PLMset objects to ExpressionSet objects. This is often useful since many Bioconductor functions operate on ExpressionSet objects. Note that the function pset2eset is a wrapper for PLMset2exprSet.

Value

```
returns a ExpressionSet
```

Author(s)

Ben Bolstad

bmb@bmbolstad.com>

See Also

ExpressionSet

Examples

```
if (require(affydata)) {
  data(Dilution)
  Pset <- fitPLM(Dilution)
  eset <- pset2eset(Pset)
}</pre>
```

preprocess

Background correct and Normalize

Description

This function pre-processes an AffyBatch.

Usage

Arguments

object an AffyBatch

subset a vector with the names of probesets to be used. If NULL then all probesets are

used.

normalize logical value. If TRUE normalize data using quantile normalization

background logical value. If TRUE background correct using RMA background correction

preprocess 13

```
background.method
name of background method to use.

normalize.method
name of normalization method to use.

background.param
list of parameters for background correction methods

normalize.param
list of parameters for normalization methods

verbosity.level
An integer specifying how much to print out. Higher values indicate more verbose. A value of 0 will print nothing
```

Details

This function carries out background correction and normalization pre-processing steps. It does not summarize to produce gene expression measures. All the same pre-processing methods supplied by threestep are supported by this function.

Value

An AffyBatch

Author(s)

Ben Bolstad

Somb@bmbolstad.com

References

Bolstad, BM (2004) Low Level Analysis of High-density Oligonucleotide Array Data: Background, Normalization and Summarization. PhD Dissertation. University of California, Berkeley.

See Also

```
expresso, rma
```

Examples

```
if (require(affydata)) {
  data(Dilution)

# should be equivalent to the bg and norm of rma()
  abatch.preprocessed <- preprocess(Dilution)
}</pre>
```

14 pseudo.coloring

· · · · · · · · · · · · · · · · · · ·

Description

These are routines used for coloring pseudo chip images.

Usage

```
pseudoPalette(low = "white", high = c("green", "red"), mid = NULL,k =50)
pseudoColorBar(x, horizontal = TRUE, col = heat.colors(50), scale = 1:length(x),k = 11, log.ticks=FA
```

Arguments

low	color at low end of scale
high	color at high end of scale
mid	color at exact middle of scale
k	number of colors to have

x A data series

horizontal If TRUE then color bar is to be draw horizontally

col colors for color bar

scale tickmarks for x if x is not numeric

log. ticks use a log type transformation to assign the colors

... additional parameters to plotting routine

Details

Adapted from similar tools in maPlots pacakge.

Author(s)

Ben Bolstad

bmb@bmbolstad.com>

See Also

```
AffyBatch, read.affybatch
```

ReadRMAExpress 15

 ${\sf ReadRMAExpress}$

Read RMAExpress computed expression values

Description

Read RMAExpress computed binary output files into a matrix or ExpressionSet

Usage

```
ReadRMAExpress(filename, return.value=c("ExpressionSet", "matrix"))
```

Arguments

filename The name of the file containing RMAExpress output to be read in

return.value should a matrix or an ExpressionSet be returned

Value

```
returns an ExpressionSet
```

Author(s)

Ben Bolstad

Somb@bmbolstad.com>

References

http://rmaexpress.bmbolstad.com

rmaPLM

Fit a RMA to Affymetrix Genechip Data as a PLMset

Description

This function converts an AffyBatch into an PLMset by fitting a multichip model. In particular we concentrate on the RMA model.

Usage

```
rmaPLM(object, subset=NULL, normalize=TRUE, background=TRUE,
    background.method="RMA.2", normalize.method="quantile",
    background.param=list(), normalize.param=list(), output.param=list(),
    model.param=list(), verbosity.level=0)
```

16 rmaPLM

Arguments

object an AffyBatch

subset a vector with the names of probesets to be used. If NULL then all probesets are

used.

normalize logical value. If TRUE normalize data using quantile normalization

background logical value. If TRUE background correct using RMA background correction

background.method

name of background method to use.

normalize.method

name of normalization method to use.

background.param

A list of parameters for background routines

normalize.param

A list of parameters for normalization routines

output.param A list of parameters controlling optional output from the routine.

model.param A list of parameters controlling model procedure

verbosity.level

An integer specifying how much to print out. Higher values indicate more ver-

bose. A value of 0 will print nothing

Details

This function fits the RMA as a Probe Level Linear models to all the probesets in an AffyBatch.

Value

An PLMset

Author(s)

Ben Bolstad <bmb@bmbolstad.com>

References

Bolstad, BM (2004) Low Level Analysis of High-density Oligonucleotide Array Data: Background, Normalization and Summarization. PhD Dissertation. University of California,

Irizarry RA, Bolstad BM, Collin F, Cope LM, Hobbs B and Speed TP (2003) *Summaries of Affymetrix GeneChip probe level data* Nucleic Acids Research 31(4):e15

Bolstad, BM, Irizarry RA, Astrand, M, and Speed, TP (2003) A Comparison of Normalization Methods for High Density Oligonucleotide Array Data Based on Bias and Variance. Bioinformatics 19(2):185-193

See Also

expresso, rma, threestep, fitPLM, threestepPLM

threestep 17

Examples

```
if (require(affydata)) {
    # A larger example testing weight image function
    data(Dilution)
    ## Not run: Pset <- rmaPLM(Dilution,output.param=list(weights=TRUE))
    ## Not run: image(Pset)
}</pre>
```

threestep

Three Step expression measures

Description

This function converts an AffyBatch into an ExpressionSet using a three step expression measure.

Usage

Arguments

object an AffyBatch.

subset a vector with the names of probesets to be used. If NULL, then all probesets are

used.

normalize logical value. If TRUE normalize data using quantile normalization

background logical value. If TRUE background correct using RMA background correction

background.method

name of background method to use.

normalize.method

name of normalization method to use.

summary.method name of summary method to use.

background.param

list of parameters for background correction methods.

normalize.param

list of parameters for normalization methods.

summary.param list of parameters for summary methods.

verbosity.level

An integer specifying how much to print out. Higher values indicate more verbose. A value of 0 will print nothing.

18 threestep

Details

This function computes the expression measure using threestep methods. Greater details can be found in a vignette.

Value

```
An ExpressionSet
```

Author(s)

Ben Bolstad <bmb@bmbolstad.com>

References

Bolstad, BM (2004) Low Level Analysis of High-density Oligonucleotide Array Data: Background, Normalization and Summarization. PhD Dissertation. University of California, Berkeley.

See Also

```
expresso, rma
```

Examples

```
if (require(affydata)) {
 data(Dilution)
 # should be equivalent to rma()
 eset <- threestep(Dilution)</pre>
 # Using Tukey Biweight summarization
 eset <- threestep(Dilution, summary.method="tukey.biweight")</pre>
 # Using Average Log2 summarization
 eset <- threestep(Dilution, summary.method="average.log")</pre>
 # Using IdealMismatch background and Tukey Biweight and no normalization.
 eset <- threestep(Dilution, normalize=FALSE,background.method="IdealMM",</pre>
                     summary.method="tukey.biweight")
 # Using average.log summarization and no background or normalization.
 eset <- threestep(Dilution, background=FALSE, normalize=FALSE,</pre>
                    background.method="IdealMM", summary.method="tukey.biweight")
 # Use threestep methodology with the rlm model fit
 eset <- threestep(Dilution, summary.method="rlm")</pre>
 # Use threestep methodology with the log of the average
 eset <- threestep(Dilution, summary.method="log.average")</pre>
 # Use threestep methodology with log 2nd largest method
 eset <- threestep(Dilution, summary.method="log.2nd.largest")</pre>
```

threestepPLM 19

```
eset <- threestep(Dilution, background.method="LESN2")
}</pre>
```

threestepPLM

Three Step expression measures returned as a PLMset

Description

This function converts an AffyBatch into an PLMset using a three step expression measure.

Usage

Arguments

object an AffyBatch

subset a vector with the names of probesets to be used. If NULL then all probesets are

used.

normalize logical value. If TRUE normalize data using quantile normalization

background logical value. If TRUE background correct using RMA background correction

background.method

name of background method to use.

normalize.method

name of normalization method to use.

summary.method name of summary method to use.

background.param

list of parameters for background correction methods

normalize.param

list of parameters for normalization methods

output.param list of parameters for output methods model.param list of parameters for model methods

verbosity.level

An integer specifying how much to print out. Higher values indicate more ver-

bose. A value of 0 will print nothing

Details

This function computes the expression measure using threestep methods. It returns a PLMset. The most important difference is that the PLMset allows you to access the residuals which the threestep function does not do.

20 threestepPLM

Value

An PLMset

Author(s)

Ben Bolstad <bmb@bmbolstad.com>

References

Bolstad, BM (2004) Low Level Analysis of High-density Oligonucleotide Array Data: Background, Normalization and Summarization. PhD Dissertation. University of California, Berkeley.

See Also

```
expresso, rma, threestep, rmaPLM, fitPLM
```

Examples

```
if (require(affydata)) {
   data(Dilution)

# should be equivalent to rma()
   ## Not run: eset <- threestepPLM(Dilution)
}</pre>
```

Index

*Topic classes	expresso, 4, 13, 16, 18, 20
PLMset-class, 9	
*Topic hplot	fitPLM, $3, 9, 16, 20$
MAplot, 5	LOUGE C. DIM. A. A. L.
*Topic manip	getCdfInfo,PLMset-method
bg.correct.LESN,2	(PLMset-class), 9
fitPLM, 3 normalize.ExpressionSet, 6 normalize.quantiles.probeset, 7 normalize.scaling, 8 PLMset2exprSet, 11 preprocess, 12 pseudo.coloring, 14	<pre>image,PLMset-method (PLMset-class), 9 indexProbes,PLMset,character-method</pre>
ReadRMAExpress, 15 rmaPLM, 15	loess, 5, 6
threestep, 17	MAplot, 5
threestepPLM, 19	MAplot, PLMset-method (MAplot), 5
AffyBatch, 2-4, 7, 8, 12-17, 19 annotation, PLMset-method (PLMset-class), 9 bg.correct.LESN, 2 boxplot, PLMset-method (PLMset-class), 9	<pre>matrix, 15 Mbox (PLMset-class), 9 Mbox, PLMset-method (PLMset-class), 9 model.description (PLMset-class), 9 model.description, PLMset-method</pre>
boxprot, refiset method (refiset crass),	mva.pairs,5
cdfName, PLMset-method (PLMset-class), 9 coefs (PLMset-class), 9 coefs, PLMset-method (PLMset-class), 9 coefs.const (PLMset-class), 9	<pre>normalize.AffyBatch.quantiles.probeset (normalize.quantiles.probeset), 7</pre>
coefs.const,PLMset-method (PLMset-class),9	normalize.AffyBatch.scaling (normalize.scaling), 8
coefs.probe(PLMset-class), 9	normalize.ExpressionSet, 6
<pre>coefs.probe,PLMset-method (PLMset-class),9</pre>	<pre>normalize.quantiles, 8 normalize.quantiles.probeset, 7 normalize.scaling, 8</pre>
<pre>coefs<- (PLMset-class), 9 coefs<-,PLMset-method (PLMset-class), 9</pre>	normvec (PLMset-class), 9 normvec, PLMset-method (PLMset-class), 9
description,PLMset-method (PLMset-class),9	NUSE (PLMset-class), 9 nuse (PLMset-class), 9
ExpressionSet, 6, 7, 9, 11, 12, 15, 17, 18	<pre>NUSE, PLMset-method (PLMset-class), 9 nuse, PLMset-method (PLMset-class), 9</pre>

22 INDEX

pData,PLMset-method(PLMset-class),9	se<- (PLMset-class), 9
pData<-,PLMset,data.frame-method	<pre>se<-,PLMset-method (PLMset-class), 9</pre>
(PLMset-class), 9	<pre>show,PLMset-method(PLMset-class),9</pre>
phenoData, PLMset-method (PLMset-class),	<pre>summary,PLMset-method(PLMset-class),9</pre>
9	, , , , , , , , , , , , , , , , , , , ,
phenoData<-,PLMset,AnnotatedDataFrame-method	threestep, 4, 13, 16, 17, 19, 20
(PLMset-class), 9	threestepPLM, 16, 19
PLMset, 3, 4, 11, 15, 16, 19, 20	,
PLMset (PLMset-class), 9	varcov (PLMset-class), 9
PLMset-class, 9	<pre>varcov,PLMset-method(PLMset-class),9</pre>
PLMset2exprSet, 11	weights (PLMset-class), 9
preprocess, 12	<pre>weights, PLMset-method (PLMset-class), 9</pre>
pset2eset (PLMset2exprSet), 11	weights<- (PLMset-class), 9
pseudo.coloring,14	<pre>weights<-,PLMset-method (PLMset-class),</pre>
pseudoColorBar(pseudo.coloring), 14	9
pseudoPalette(pseudo.coloring),14	
quantile, 8	
read.affybatch, 14	
ReadRMAExpress, 15	
resid(PLMset-class),9	
resid,PLMset-method(PLMset-class),9	
resid<-(PLMset-class),9	
resid<-,PLMset-method(PLMset-class),9	
residSE(PLMset-class),9	
residSE,PLMset-method(PLMset-class),9	
residuals (PLMset-class), 9	
residuals, PLMset-method (PLMset-class),	
9	
residuals<- (PLMset-class), 9	
residuals<-,PLMset-method	
(PLMset-class), 9	
RLE (PLMset-class), 9	
RLE, PLMset-method (PLMset-class), 9	
rma, 4, 13, 16, 18, 20	
rmaPLM, 15, 20	
compleNemes DI Most method	
sampleNames,PLMset-method (PLMset-class),9	
, , , , , , , , , , , , , , , , , , , ,	
sampleNames<- (PLMset-class), 9	
sampleNames<-,PLMset,character-method	
(PLMset-class), 9	
se (PLMset-class), 9	
se, PLMset-method (PLMset-class), 9	
se.const (PLMset-class), 9	
se.const,PLMset-method(PLMset-class),9	
se.probe(PLMset-class),9	
se.probe,PLMset-method(PLMset-class),9	